**X-RAY TUBE
HOUSING ASSEMBLY**

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RÖNTGENSTRAHLER

***EMISOR DE
RAYOS X***

SRO 2550 ROT 350

**INSTRUCTIONS FOR USE / SERVICE INFORMATION
INSTRUCTIONS D'UTILISATION / INFORMATION DE SERVICE
GEBRAUCHSANWEISUNG / SERVICEINFORMATION
*INSTRUCCIONES DE USO / INFORMACIÓN DE SERVICIO***

**English
Français
Deutsch
*Español***

PHILIPS

Copyright

© 2004 Philips Medical Systems
DMC GmbH

All rights reserved. Reproduction in whole or in part is prohibited without the prior written consent of the copyright-holder.

Philips Medical Systems reserves the right to make changes in specifications or to discontinue any product, at any time without notice or obligation, and is not liable for any consequences resulting from the use of this publication.

4512-984-28101 REV AC - 03/2005

Internet address: www.medical.philips.com

Equipment is subject to change without notice. All changes will be in compliance with regulations governing manufacture of medical equipment. Printed in Germany.

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INSTRUCTIONS FOR USE / SERVICE INFORMATION

X-ray tube housing assembly

SRO 2550 ROT 350**9890-000-85831****SRO 2550 ROT 350 180 °****9874-004-24112**

In case there are any questions concerning this manual,
please send this LOPAD via fax to (+49) 40/50 78-24 81.

File: SRO_2550_ROT_350_28101_AC

List of pages and drawings (LOPAD)

Manual Order No.: 4512-984-28101
released: 03/2005

01 ... 06		Z-3	(04.0)	A4	SRO 2550
<hr/>					
1 ... 28	(05.0)	English	Z-4	(04.0)	A4 SRO 2550
<hr/>					
1 ... 28	(05.0)	French	Z-5.1	(04.0)	A4 SRO 2550
<hr/>					
1 ... 28	(05.0)	French	Z-5.2	(04.0)	A4 SRO 2550
<hr/>					
1 ... 28	(05.0)	German	Z-6	(04.0)	A4 SRO 2550
<hr/>					
1 ... 28	(05.0)	German	Z-7.1	(04.0)	A4 SRO 2550
<hr/>					
1 ... 28	(05.0)	Spanish	Z-7.2	(04.0)	A4 SRO 2550
<hr/>					
		Z-8	(04.0)	A4	SRO 2550
Z-1.1	(04.1)	A4	Z-10	(04.0)	A4 ROT 350 ... ROT 360
Z-1.2	(04.1)	A4	Z-11	(04.1)	A4 ROT 350 Ali Ring
Z-1.3	(04.1)	A4	Z-12	(04.1)	A4 ROT 350; ROT 360
Z-1.4	(04.1)	A4	Z-13	(04.0)	A4 ROT 350; ROT 360
Z-2.1	(05.0)	A4	<hr/>		
Z-2.2	(05.0)	A4			
Z-2.3	(05.0)	A4			
Z-2.4	(05.0)	A4			
			<hr/>		
			OTR 4512-988-0003x		

INSTRUCTIONS FOR USE / SERVICE INFORMATION

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11

12

INSTRUCTIONS FOR USE / SERVICE INFORMATION**TEXT**

	CONTENTS	1
1	BRIEF DESCRIPTION OF THE X-RAY TUBE HOUSING ASSEMBLY.....	3
2	SAFETY INFORMATION	5
2.1	Radiation protection	6
2.2	Electrical safety	6
2.3	Thermal safety	6
2.4	Electromagnetic compatibility	7
2.5	Disposal	7
2.6	Conformity	8
3	COMPATIBILITY.....	9
3.1	Compatible generators	9
3.2	Compatible rotor control units	9
4	TECHNICAL INFORMATION	10
4.1	Denotation	10
4.2	Identification	10
4.3	Accessories	10
4.4	Technical data	11
4.5	Filtration	11
4.6	Collimation near the focus	11
4.7	Position of the cable sockets with respect to the central beam	11
4.8	Universal spacing plates	12
4.9	Temperature switch	12
4.10	HEW lamp	12
4.11	Fan	12
5	INSTALLATION	13
5.1	Mechanical installation	13
5.1.1	Alignment ring	13
5.1.2	Trunnion ring	14
5.2	Electrical connections	15
5.2.1	Connection of the protective-earth conductor	15
5.2.2	Connection of the high voltage cables	15
5.2.3	Connection of the stator cable	17
5.2.4	Connection of the fan	18
5.2.5	Connection of the excess temperature switch	18
6	SETTING-TO-WORK	19
6.1	X-ray tube conditioning	19
6.2	Additional information	21
7	MAINTENANCE	22
7.1	Corrective maintenance	22
7.2	Planned maintenance	22
7.2.1	Maintenance by the user	25
7.2.2	Maintenance by a service organization	26
7.3	Return shipment process	27

TECHNICAL DATA

Technical data of X-ray tube	Z-1.1
Technical data of X-ray tube housing assembly	Z-2.1
Emission and filament curves	Z-3
Charts for single load rating	Z-4
Heating and cooling curves: Anode of X-ray tube	Z-5.1
Heating and cooling curves: X-ray tube housing assembly	Z-5.2
Load data: Cine radiography	Z-6
Load data: Serial radiography	Z-7.1
Load data: Serial radiography	Z-7.2
Mechanical dimensions of X-ray tube	Z-8
Labelling of X-ray tube housing assembly	Z-10
Mechanical dimensions of X-ray tube housing assembly	Z-11
Electrical connections	Z-12
Installation of fan	Z-13

1 BRIEF DESCRIPTION OF THE X-RAY TUBE HOUSING ASSEMBLY

The high performance SRO 2550 rotating anode X-ray tube is integrated into the X-ray tube housing assembly.

X-ray tube housing assembly

See
figure 1

1. Radiation protection housing
2. X-ray tube
3. Connector plate for the collimator
4. X-ray tube window with filter zone


Figure 1: X-ray tube housing assembly

X-ray tube

See
figure 2

5. Cathode
6. Anode disk
7. Rotor


Figure 2: X-ray tube

Due to the high energy dissipation of the anode and to the large X-ray field, this X-ray tube is suitable for a universal field of application.

It has mainly been developed for applications in radiography:

- Serial exposure
- Fluoroscopy technique

The SRO 2550 X-ray tube is installed in the X-ray tube housing assembly:

SRO 2550 ROT 350

SRO 2550 ROT 350 180 °

The unit is only available as a complete X-ray tube housing assembly.

The exchange of an X-ray tube is possible only at the Hamburg factory and at the reloading stations of Philips Medical Systems.

2 SAFETY INFORMATION

This manual is designed to make it possible for you to work safely with the X-ray tube housing assembly described. You may operate the X-ray tube housing assembly only in compliance with the safety instructions in this manual and not use it for purposes other than for which it is intended. The X-ray unit may only be operated by persons who have the necessary expertise in radiation protection and who have been instructed how to operate the X-ray unit.

It is always the user who is responsible for compliance with the regulations applying to installation and operation of X-ray units.

- You must never use the X-ray tube housing assembly if it has any electrical, mechanical or radiological defects. This particularly applies to faulty indicators, displays, warnings and alarms.
- If the user wishes to connect the X-ray tube housing assembly to other equipment, components or assemblies and if it is not apparent from the technical data whether it can be safely combined with such equipment, components or assemblies, the user must ensure that the safety of the
 - patient
 - operating staff
 - third parties
 - environment

is not affected by the planned combination by consulting the manufacturers involved or by making enquiries from an expert.

- Philips is responsible for the safety features of its products only if maintenance, repairs and modifications have been performed by the service organization or by persons explicitly authorized to do so by the manufacturer.
- As with any technical appliance, this equipment requires
 - correct operation
 - regular, competent maintenance
 - care

which are described in chapter 7 "MAINTENANCE".

- If you operate the respective X-ray equipment, and hence the X-ray tube housing assembly incorrectly or if the user fails to maintain it properly, the manufacturer cannot be held responsible for any resulting malfunctions, damage or injuries.
- The safety protection switch for the X-ray tube housing assembly, which prevents radiation from being switched ON if maximum X-ray tube housing assembly temperature is exceeded, must be neither removed nor bypassed.
- To ensure patient safety, some generators allow emergency fluoroscopy even though the X-ray tube housing assembly excess temperature switch has been activated.
- You may remove or open parts of the housing only if you are instructed to do so in this manual.

2.1 Radiation protection

RADIATION

Make certain before every X-ray exposure that all necessary radiation precautions have been taken.

You can find information about radiation precautions in the instructions for use for the individual X-ray systems with which you are using this X-ray tube housing assembly.

2.2 Electrical safety

WARNING

Only trained maintenance staff may remove the covers from the X-ray tube housing assembly.

This X-ray tube housing assembly may only be used in medical rooms which meet the requirements of relevant national and international standards and laws.

- You must not use this X-ray tube housing assembly in areas where there is danger of explosion.
- Detergents and disinfectants, including those used on the patient, can create explosive gas mixtures.

Please observe the relevant regulations.

2.3 Thermal safety

WARNING

To prevent damage due to thermal overloading, precautions must be taken to ensure that the X-ray tube housing assembly is not operated outside its specified load parameters. This is the only way to eliminate all risk to the patient and operator.

The X-ray tube housing assembly, consisting of the X-ray tube and X-ray tube housing, is one component in the X-ray system. The power supply and the frame rate are generated and controlled by the X-ray system. The user interfaces of some X-ray generators indicate the thermal state of the X-ray tube.

2.4 Electromagnetic compatibility

EMC

In accordance with its intended use, this electronic product complies with the law governing EMC, which defines the permitted emission levels from electronic equipment and its required immunity against electromagnetic fields.

Electronic apparatus that satisfies the EMC requirements is designed so that under normal conditions there is no risk of malfunction caused by electromagnetic interference. However, in the case of radio signals from high-frequency transmitters with a relatively high transmitting power, the risk of electromagnetic incompatibility when operated in close proximity to electronic apparatus cannot be totally ruled out.

WARNING

In unusual circumstances unintended functions of the apparatus could be initiated, possibly giving rise to undesirable risks for the patient or user. For this reason, all kinds of transmission with mobile radio equipment should be avoided. This also applies when the apparatus is in STANDBY mode. Mobile telephones must be switched OFF in designated problem zones.

2.5 Disposal

Philips manufactures state-of-the-art X-ray equipment in terms of safety and environmental protection. Assuming no parts of the system housing are opened and assuming the system is used properly there are no risks to persons or the environment.

To comply with regulations it is necessary to use materials which may be harmful to the environment and therefore have to be disposed of in a proper manner.

WARNING

This X-ray tube housing assembly contains materials which are toxic. For this reason you must not dispose of the X-ray equipment together with industrial or domestic waste.

Philips

- supports you in disposing of the X-ray equipment described in a proper manner.
- recycles used materials and components from X-ray tubes and X-ray tube housing assemblies. Elaborate testing and quality assurance procedures, including detailed inspection of sub-components, ensures the same high levels of quality and functionality provided by new materials.
- helps to reduce environmental pollution.

Consequently, do contact your Philips Sales Organization in full confidence.

2.6 Conformity

CE

The Philips X-ray tube housing assembly described in this manual meets the provisions of the Medical Device Directive 93/42 EEC (93).

This manual was drafted in Germany.

If you have further questions regarding legal issues and applicable national or international standards, please address them to:

Philips
Medical Systems DMC GmbH
Röntgenstr. 24
D-22335 Hamburg

Fax: (+49) 40/50 78-21 47

3 COMPATIBILITY

The X-ray tube housing assembly is compatible with:

- Generators
- Rotor control units
- Collimators

3.1 Compatible generators

- | | |
|---------------|---------------|
| – MEDIO CP | 30 / 50 |
| – MEDIO CP-H | 50 / 65 |
| – SUPER CP-D | 50 / 80 |
| – SUPER CP | 50 / 80 / 100 |
| – OPTIMUS CP | 1050 / 2000 |
| – OPTIMUS RAD | 50 / 65 / 80 |
| – OPTIMUS R/F | 50 / 65 / 80 |
| – OPTIMUS C | 50 / 65 / 80 |
| – Velara R/F | 65 / 80 / 100 |

3.2 Compatible rotor control units

- MEDIO SRO rotor control for D 88
- MEDIO SRO rotor control - H - *
- MEDIO SRO rotor control *
- **RO** (/ SRO / SRM) modular rotor control **50 / 60 Hz** (150 / 180 Hz) **
- (RO /) **SRO / SRM** modular rotor control (50 / 60 Hz) **150 / 180 Hz**
- SRC rotor control
- SCP-RO rotor control **
- SCP SRO + SRM rotor control module 517
- Rotor control dual speed rotation

*) Without fluoroscopy rotation, without fast exposure.

**) Only for fluoroscopy rotation.

4 TECHNICAL INFORMATION

4.1 Denotation

The name of the X-ray tube housing assembly means:

SRO	Super Rotalix X-ray tube, high speed rotation
25	Nominal anode input power for small focal spot [kW]
50	Nominal anode input power for large focal spot [kW]
ROT 350	X-ray tube housing of Philips family 350, air cooled, 90 ° (position of the cable sockets with respect to the central beam)
ROT 350 180 °	X-ray tube housing of Philips family 350, air cooled, 180 ° (position of the cable sockets with respect to the central beam)

4.2 Identification

Type numbers X-ray tube housing assembly

SRO 2550 ROT 350	9890-000-85831
SRO 2550 ROT 350 180 °	9874-004-24112

Type numbers X-ray tube

SRO 2550	9806-206-11102
----------	----------------

4.3 Accessories

– HEW lamp	Type No.: 4512-148-2552x
– Fan	Type No.: 9806-705-0020x

4.4 Technical data

For technical data of the X-ray tube refer to drawing Z-1.

For technical data of the X-ray tube housing assembly refer to drawing Z-2.

4.5 Filtration

The minimum total filtration value of the X-ray tube housing assembly is 2.5 mm Al equivalent (related to half the rated voltage according to IEC 60522).

Thus it complies with international regulations.

An additional filter plate of 1.7 mm Al which can only be removed by special tools is part of total filtration.

Inherent filter value: 0.8 mm Al equivalent.

WARNING

It is only allowed to remove the additional filter plate of 1.7 mm Al when the manufacturer of the used X-ray system guarantees a minimum total filtration of 2.5 mm Al equivalent.

The filter value is indicated on the type label of the X-ray tube housing assembly.

4.6 Collimation near the focus

Collimation near the focus is effected by a lead aperture plate which depends upon the X-ray tube installed since the anode angles of the different X-ray tubes may vary. This lead plate reduces scattered radiation.

Caution!

In case the field of radiation actually needed for applications is smaller than the one limited by the lead plate installed by the factory, the lead plate has to be exchanged. Philips has to be contacted.

4.7 Position of the cable sockets with respect to the central beam

The position of the cable sockets with respect to the central beam is defined with 90 ° (standard value). If it is not 90 °, the position is defined by an abbreviated designation, e.g. "180 °".

It reflects the angle between cable socket and central beam, measured clockwise, with the X-ray tube housing viewed from the cathode side and the cable sockets pointing vertically upwards.

4.8 Universal spacing plates

These plates are located at the X-ray port and the mounting surface of the collimator. Some intermediate spacing plates are placed by the factory depending on the respective focus used.

Caution! The factory adjustment must not be changed during installation work.

4.9 Temperature switch

The temperature switch is installed at the cathode side of the X-ray tube housing assembly. It is activated by oil temperatures in the X-ray tube housing at about 85 °C. This temperature switch interrupts the READY state of the system.

4.10 HEW lamp

The lamp is part of accessories. It does not belong to the X-ray tube housing assembly.

This lamp indicates which X-ray tube housing assembly is selected and whether it is ready for operation.

If required by the respective national regulations, this lamp has to be installed by the local field service.

4.11 Fan

The fan is part of accessories. It does not belong to the X-ray tube housing assembly.

With fan installed, the max. continuous heat dissipation of the X-ray tube housing assembly can be increased to 350 W.

The fan increases the weight of the X-ray tube housing assembly by 1 kg.

5 INSTALLATION

5.1 Mechanical installation

The attachment of the X-ray tube housing assembly is described in the installation instructions of the system.

5.1.1 Alignment ring

See figure 3 The example in the figure beside shows the SRO-family.

The four screws of the alignment ring are sealed to guarantee the correct fit.

The collimator is automatically aligned in the correct position.


Figure 3: Sealing of the alignment ring

If a compatible collimator is used

Caution! The alignment ring must not be removed.

- Mount the coupling flange on the alignment ring respective the correct position of the fit, slot and key.
- The collimator can be mounted without further alignment.

If a not compatible collimator is used

The alignment ring is mounted when the tube is delivered.

- Remove the alignment ring and dispose it.
It is not recyclable.
- Proceed the installation in the common way.

See
figure 4

1. Coupling flange
2. Alignment ring


Figure 4: Removal of the alignment ring

5.1.2 Trunnion ring

See
drawing Z-11

The trunnion ring of the X-ray tube housing holds the X-ray tube housing assembly in two clamp rings. It provides that the X-ray tube housing assembly can rotate around its longitudinal axis. This movement can be blocked by tightening the two clamp rings.

The adjustment of the central beam in longitudinal and transverse direction with respect to the fixing of the X-ray tube housing assembly must be carried out on the system.

5.2 Electrical connections

5.2.1 Connection of the protective-earth conductor

Note: Use the special material in the accessories supplied.

The X-ray tube housing assembly must be earthed with a separately routed protective-earth (PE) conductor cable of at least 4 mm² Cu.

See
drawing Z-12

- Connect the PE conductor to the central earthing point of the X-ray tube housing assembly. Use the ring terminal in the accessories supplied.
- Establish connection between the central earthing point of the X-ray tube housing assembly and the cover cap of the X-ray tube housing.

5.2.2 Connection of the high voltage cables

WARNING

Discharge the H.V. cables after having pulled them out. They act like a capacitor.

Note: Use old plug receptacles R3 only in combination with adapter O3/R3, type No. 9806-420-4000x.
In such case the maximum length of each high voltage cable is 20 m.
For cine operation special restrictions are applicable.

See
figure 5

The connector pins of the high voltage cables are slotted and the two legs of the pins must have specified distance with respect to each other. If this specified distance is changed, i.e. if it is reduced, the legs of the pins do not longer touch the inner surface of the contact jacks such that contact is ensured. In this case, minor voltage sparkover may occur which burns or impairs the contact surface.

It might happen during installation work or an exchange of the X-ray tube housing assembly that the connector is exposed to too much mechanical stress and that the legs of the contact are pressed together.

- In this case readjust gap of splint pins. Use pin gage tool, type No. 9890-000-0284x.


Figure 5: Conditions of the connector pins

Caution! When connecting the high voltage cables to the X-ray tube housing assembly, use the silicone washers supplied and use nothing but silicone paste, even when exchanging the X-ray tube housing assembly.

- Provide the silicone washer with a thin layer of silicone paste.
- Fit the silicone washer over the pins of the high voltage plug. The connector pins must be free from silicone paste. When connecting the high voltage cables to the high voltage generator, do not use the silicone washers and the rubber ring.
- Check whether the high voltage plug is clean.
- If necessary, clean it with alcohol. Do not use anything else like benzine or trichlorethylene.
- Fit the high voltage plug. While doing so, observe the position of the nipple.
- Tighten the coupling ring and secure it.

5.2.3 Connection of the stator cable

Use a stator cable 7 x 0.82 mm².

Another type of stator cable 3 x 1.31 mm², screened has to be used only for connecting the stator to the rotor control unit of a generator of the OPTIMUS/Velara family.

Note: *If using a generator of the OPTIMUS/Velara family earth the screening of the stator cable with a metallic clamp.*

Rotor control SRO


Figure 6: Stator connection

*) Generator of the OPTIMUS/Velara family only.

- See figure 6
- Place the jumpers across terminals XA 100 and XA 200.
 - Connect the stator cable:

Rotor control	Cable wire	Phase	Connection terminal	Stator data phase/phase
SRO	1	U	XA 101	U - V: 11 Ω ± 10% 57 mH ± 10% V - W: 9 Ω ± 10% 34 mH ± 10%
	2	V	XA 109	
	3	W	XA 202	

Table 1: Stator connection

5.2.4 Connection of the fan

Stator cable 7 x 0.82 mm²

Use two wires of the stator cable.

3-phase stator cable

An additional cable 2 x 0.82 mm², screened has to be installed at the feedthrough of the stator cable.

See
drawing Z-13

- Remove the protective cap from the anode side.
- Fix the bracket (2) with two countersunk screws M3 to ring (3).
- Arrange the leads going to the stator and fan connection and squeeze them down so far that guard (4) can be placed upon ring (3). Fix the guard to the ring with two cheese-head screws M3, washers and lock washers.
- Connect the leads to the fan (2 x 0.82 mm²).
The operating voltage 230 V of the fan can be taken from the 230 V supply circuit of the generator.
- Place cap (5) with the built-in fan upon the X-ray tube housing in such a manner that the three through-holes at the periphery of the cap coincide with the three threaded holes M5 on ring (3). One of the three threaded holes is located in bracket (2). Fix the cap with three special-type screws (1) and three locking washers.
- Before attaching the anode and cathode end caps mechanically, connect them with the protective-earth conductor.
- Bring end cap (6) into position and fix it with three countersunk screws.
- Check that the vane of the fan is not obstructed.

5.2.5 Connection of the excess temperature switch

Note: Use a 2 x 0.82 mm² screened cable.

See
drawing Z-12

- Solder the cable on at the excess temperature switch.
- Secure the cable with a cable clamp at the X-ray tube housing to provide drag relief.

6 SETTING-TO-WORK

Note: *Modular generators of the small and big family made in Hamburg are out of service. For information about installation of an X-ray tube housing assembly together with modular generator it is necessary to contact the Helpdesk X-ray in Hamburg.*

Setting-to-work includes:

- Software and hardware programming for the X-ray tube type.
- Adjustment for the filament.
- Adaptation and running-in of the X-ray tube.

6.1 X-ray tube conditioning

RADIATION

**Radiation is released during the conditioning process.
Observe the regulations for radiation protection.**

The conditioning process must be performed for each new X-ray tube housing assembly to be connected, irrespective of the storage time.

The intervals between exposures must be adhered too and monitored by a watch.

Note: *Conditioning has to be performed with the **large focal spot** only.*

The first exposure of each day does not require warming up.

For operating breaks of more than 3 months, an exposure-series have to be made with the specified kV and mAs values as given in the table below.

The generator must be in READY state, i.e. the green LED on the desk must be lit. Pay attention that the READY state conditions are met.

The exposure run (kV - mAs / kV - mAs - s and break time) must always be performed at the selected kV level without repeating the start-up procedure of the X-ray tube housing assembly, i.e. the PREP button must remain to be pressed during the run.

Conditioning procedure

- Switch ON the generator and select an appropriately programmed auxiliary for the X-ray tube housing assembly to be conditioned.
- Select the large focal spot of the X-ray tube housing assembly.
- Select 2-button-technique "kV-mAs".
Choose appropriate values at the generator desk corresponding to the following table.
- Release exposures according to the following table.

Note: *In case of electrical interferences, the procedure must be continued after an interval of 5 min, recommencing at a lower kV level.*

Number of exposures	X-ray tube voltage [kV] *	mAs product [mAs]	Break time [s]
5	81	125	1
			30
3	102	125	1
			30
2	117	125	1
			60
2	141	125	1
			120

Table 2: Values for conditioning of the X-ray tube housing assembly

*) Given values depend on the generator.
They may be slightly different.

6.2 Additional information

Measurement of X-ray tube high voltage and X-ray tube current

Use only the appropriate instruments for measuring the X-ray tube high voltage and X-ray tube current.

Exposure series

The tables Z-7 give the maximum number of exposures permissible per series.
The minimum cooling time between subsequent series is 20 min.

Cine radiographs

The tables Z-6 give the max. anode input power of each radiography in kW based on equivalent anode input power of 250 W.

7 MAINTENANCE

7.1 Corrective maintenance

Repairs of X-ray tube housing assemblies must be carried out only in the X-ray tube factory.

Exception: The provider has a properly equipped service workshop.

After installation of a new X-ray tube housing assembly:

- Check the high voltage cables and use new silicone washers.
See chapter 5.2.2 "Connection of the high voltage cable".
- Exchange the identification label of the X-ray tube housing assembly at the label bracket of the system.

7.2 Planned maintenance

As with any technical appliance these X-ray tube housing assemblies also require:

- Regular testing by the user.
- Regular service and repair.

By taking these precautions you maintain the operability and operational reliability of the system. As the user of an X-ray unit you are obliged according to accident standard regulations, the medical products law and other regulations to take such precautions.

Maintenance consists of tests which the user can perform and maintenance which is performed by qualified and trained specialists.

WARNING

X-ray tube housing assemblies contain mechanical components which are subjected to wear and tear due to operation.

The correct setting of the electromechanical and electronic assemblies safeguards the functioning, image quality, electrical safety and exposure of the patient and medical personnel to radiation.

Philips recommends you to:

- Perform the tests indicated in the tables of chapter 7.2.1 "Maintenance by the user" on a regular basis.
- Have the X-ray unit serviced by a service organization at least once a year.
You must have heavily used X-ray equipment subjected to maintenance more frequently.

In this way you avoid endangering the patient and you meet your obligations.

By entering into a service agreement with a service organization you retain the value and safety of your X-ray equipment. All the necessary maintenance, including the safety tests for the purpose of standard avoidance of danger and the necessary settings for optimum image quality and minimum exposure to radiation, are performed at regular intervals. The service organization agrees on these intervals with you, taking the legal requirements into account.

If operational defects or other departures from normal operational behavior occur, you must switch OFF the X-ray unit and inform your respective service organization. You may only resume operation of the X-ray equipment when it has been repaired. Operation using faulty components may lead to an increased safety risk or unnecessarily high exposure to radiation.

WARNING

Faulty components which affect the safety of the X-ray equipment must be replaced by genuine spare parts.

Recording results

Service and repairs must be entered in the medical products logbook, including the following data:

- Type and scope of work.
- Details of any changes to ratings or the working zone.
- Date, company and person performing the work, signature.

Cleaning

Detergents with a high alcohol content make the material dull or cause it to crack.

Caution!

**Switch OFF the mains voltage of the X-ray equipment before cleaning.
Never use any caustic, solvent or abrasive detergents or polishes.**

Heed during cleaning:

- Be careful that no water or other fluids can get inside the X-ray equipment.
This prevents short circuits in the electrical installations and corrosion on the components.
- Clean enamelled parts and aluminium surfaces only with a moist cloth and a mild detergent and wipe with a dry woollen cloth.
- Wipe chrome parts only with a dry woollen cloth.

Disinfection The disinfection method used must conform to valid legal requirements and guidelines for disinfection and protection against explosion.

WARNING

Switch OFF the mains voltage before disinfection of the X-ray equipment. If you use disinfectants which create explosive gas mixtures, these must first have dissipated before restarting the X-ray equipment.

Wipe disinfection

Disinfect all parts of the X-ray equipment, including accessories and connecting cables, by wiping only with a soft cloth.

Spray disinfection

Spray disinfection is not recommended because disinfectant could get inside the X-ray equipment.

Disinfection by atomizing

Cover the X-ray equipment carefully with a sheeting when it has cooled down. After the disinfection mist has subsided, you may remove the sheeting and disinfect the X-ray equipment by wiping with a soft cloth.

7.2.1 Maintenance by the user

- Visual check**
- Check the X-ray equipment for apparent defects:

Interval	Scope of work	Method
As per relevant national and international standards and laws or the local regulations.	Stability test	As per relevant national and international standards and laws or the local regulations.
Daily	Damaged parts, labelling and warning plates	Inspection
Weekly	All cables and connections (damage, breaks)	Inspection
Weekly	Oil leaks and unusual noises	Inspection

Table 3: Defect check

- Error messages**
- Check the X-ray equipment for error messages:

Type	Meaning	Action
Warning	The X-ray tube housing assembly is thermally overloaded and cannot operate further without a rest period.	Maintain rest intervals and wait for enabling at generator desk.

Table 4: Error message check

In case of other error messages concerning the X-ray tube housing assembly call your local service.

7.2.2 Maintenance by a service organization

The following maintenance work must be carried out **once a year**.

- Check the X-ray tube housing assembly for obvious damage that might impair radiation protection and may affect safety.
Example: Dents, scratches, tears, soiling, wear, oil leaks, unusual noises, ...
- Check and confirm that the Al filter provided is properly located and fixed in the useful beam.
- Check all screws of the X-ray tube housing assembly which determine the position.
- Check all cables and connections for mechanical damage, burned areas, loosened connectors and the completeness of electrical connections.
Example: Protective earth, high voltage cables, temperature switch, stator cable, ...
- Check the necessary monitoring, safety, display and indicating systems.
If two or more X-ray tube housing assemblies are controlled by a single exposure switch, check whether the working site selected is indicated by the lamp mounted on or in the vicinity of the X-ray tube housing assembly.
Observe the local national regulations.
Example: X-ray ON, selected X-ray tube, X-ray tube load indication, ...

Safety checks regarding national and international standards and laws or the local regulations

- Measuring the safety-relevant output parameters.
Example: Stability test, ...
- For the particular product other special technical tests according to the generally accepted standards of engineering practice.

Data logging

- Record results and test reports in the X-ray system manual (Medical Products Logbook).

7.3 Return shipment process

Use the re-usable packing of the new delivered X-ray tube housing assembly for the return shipment of the faulty X-ray tube housing assembly.

This guarantees a transport of the defective X-ray tube housing assembly in a packing which has been developed for this purpose.

All X-ray tube housing assemblies are delivered including a set of documentation. In addition to the Instructions For Use / Service Information and the CE certification this set of documentation also contains a Unit Manual "Online Tube Report" which describes the procedure to generate the Tube Report. See section 8.

The Tube Report also consists of a CD-ROM inclusive software filed at the beginning of this manual.

After an exchange of the X-ray tube housing assembly:

- Proceed the instructions of the Unit Manual "Online Tube Report" of section 8.
- Return the defective X-ray tube housing assembly to the X-ray tube factory in Hamburg.

Address: **Philips
Medical Systems DMC GmbH
Supply & Service Logistics / Returns
Building G
Röntgenstr. 24
D-22335 Hamburg**

- Send the Tube Report via E-Mail to DMC Hamburg, BU X-Ray Tubes.

E-Mail address: **ms.de.tubes.datalog@philips.com**

When the X-ray tube housing assembly arrives in Hamburg and the Tube Report is completely filled out the customer receives a credit note for the X-ray tube housing or warranty is issued according to the warranty regulations arranged between BU X-Ray Tubes and Customer Services.

Also see Intranet: **http://pww.xray.ms.philips.com/cmt_xray/
Products ---> X-ray Tubes ---> Warranty Regulation**

Note regarding the X-ray tube housing credit note

In case the Tube Report is missing or not completely filled out the credit note for the X-ray tube housing is reduced respectively!

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
1	Nominal X-ray tube voltage Symmetrical with respect to ground, 6- and 12- pulse equivalent or DC	kV	150	
2	Anode material Rhenium alloyed tungsten compound anode: RT / TZM	---		
3	Location of focal spots: Superimposed	---		
4	Nominal focal spot values, superimposed (IEC 60336)	---	0.6	1.0
5	Nominal anode input power, based on 0.1 s (IEC 60613) 20 W equivalent anode input power 250 W equivalent anode input power	kW	30 25	60 50
6	Max. anode heat content	kJ kHU	220 300	
7	Max. continuous heat dissipation Rotating anode Non-rotating anode	W	450 350	
8	Max. anode heat dissipation (Not for continuous operation)	W HU/min	1300 105300	
9	Anode disk diameter	mm	90	
10	Anode angle	°	15	
11	Max. usable X-ray field in 100 cm SID	cm	50 x 50	
12	Nominal acceleration time	s	1.0 *	
13	Nominal braking time	s	approx. 10.0 *	
14	Nominal anode speeds	RPM	3000 / 3600 9000 / 10800	
15	Momentum of inertia of rotor with anode	g x cm ²	6150	
16	Sense of anode rotation (Seen from cathode side)	---	counter-clockwise	
17	Filament diameter	µm	250	
18	Max. filament current CP generator OPTIMUS/Velara generator	A	5.79 ** 6.26 **	6.49 ** 6.54 **
19	Max. filament voltage CP generator OPTIMUS/Velara generator	V	10.2 ** 10.9 **	15.0 ** 14.6 **
20	Weight of tube	kg	2.25	

*) Depending on rotor control unit.

**) Values are for adaptation and also depending on generator control.

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



Note: Emission characteristics valid for 6/12-pulse and converter generators.

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
Note: Based on 6-/12-pulse or DC and equivalent anode input power 250 W at high speed rotation.

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Nota: Temperatura ambiental $\leq 25\text{ }^{\circ}\text{C}$

Heating and cooling curves: Anode of X-ray tube

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

(04.0)

Z-5.1

Note: Including heat content of the anode, heat dissipation of the stator and cathode supply at ambient temperature $\leq 25^{\circ}\text{C}$.

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Heating and cooling curves: X-ray tube housing assembly

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The tables give the max. anode input power of each radiography in kW based on equivalent anode power of 250 W.

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The tables give the max. number of exposures permissible per series.
The min. cooling time between subsequent series is 20 min.

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(04.0)

Z-7.1

© 2004 Philips Medical Systems
ALL RIGHTS RESERVED

Focal spot: 1.0			RPM: 9000/10800					Anode input power: 50 kW				
exp./s	Exposure time of individual radiograph [ms]											
	3	6	10	20	32	40	50	63	80	100	120	160
12	100	100	100	58	21	11	4	-	-	-	-	-
10	100	100	100	68	27	15	8	3	-	-	-	-
8	100	100	100	81	34	21	12	6	-	-	-	-
6	100	100	100	97	44	29	18	11	5	-	-	-
4	100	100	100	100	58	40	27	17	10	6	-	-
2	100	100	100	100	85	60	43	30	20	13	9	5
1	100	100	100	100	100	84	61	43	30	21	16	9

Focal spot: 1.0			RPM: 9000/10800					Anode input power: 40 kW				
exp./s	Exposure time of individual radiograph [ms]											
	3	6	10	20	32	40	50	63	80	100	120	160
12	100	100	100	100	37	22	12	5	-	-	-	-
10	100	100	100	100	45	28	17	8	3	-	-	-
8	100	100	100	100	55	36	23	13	6	-	-	-
6	100	100	100	100	68	46	31	19	11	6	-	-
4	100	100	100	100	87	61	42	28	18	11	7	3
2	100	100	100	100	100	89	63	44	30	21	15	9
1	100	100	100	100	100	100	88	63	44	31	24	15

Focal spot: 1.0			RPM: 9000/10800					Anode input power: 30 kW				
exp./s	Exposure time of individual radiograph [ms]											
	3	6	10	20	32	40	50	63	80	100	120	160
12	100	100	100	100	72	46	28	16	7	-	-	-
10	100	100	100	100	83	55	35	21	11	5	-	-
8	100	100	100	100	97	66	44	28	16	9	5	-
6	100	100	100	100	100	81	55	37	23	14	9	3
4	100	100	100	100	100	100	72	49	33	22	15	8
2	100	100	100	100	100	100	100	73	51	36	26	16
1	100	100	100	100	100	100	100	100	71	51	39	25

The tables give the max. number of exposures permissible per series.

The min. cooling time between subsequent series is 20 min.

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